

Health & Care Information Model:

nl.zorg.LaboratoryTestResult-v4.0

Status:Final
Release:2017
Release status: Prepublished

Managed by:

Better health
through better IT



Content

| | |
|---|----------|
| 1. nl.zorg.LaboratoryTestResult-v4.0 | 3 |
| 1.1 Revision History | 3 |
| 1.2 Concept | 4 |
| 1.3 Mindmap | 4 |
| 1.4 Purpose | 4 |
| 1.5 Patient Population | 4 |
| 1.6 Evidence Base | 4 |
| 1.7 Information Model | 4 |
| 1.7.1 LaboratoriumTest | 7 |
| 1.7.2 Monster | 12 |
| 1.8 Example Instances | 17 |
| 1.9 Instructions | 17 |
| 1.10 Interpretation | 18 |
| 1.11 Care Process | 18 |
| 1.12 Example of the Instrument | 18 |
| 1.13 Constraints | 18 |
| 1.14 Issues | 18 |
| 1.15 References | 18 |
| 1.16 Functional Model | 18 |
| 1.17 Traceability to other Standards | 18 |
| 1.18 Disclaimer | 18 |
| 1.19 Terms of Use | 19 |
| 1.20 Copyrights | 19 |

1. nl.zorg.LaboratoryTestResult-v4.0

| | |
|---------------------------------|--|
| DCM::CoderList | Kerngroep Registratie aan de Bron |
| DCM::ContactInformation.Address | * |
| DCM::ContactInformation.Name | * |
| DCM::ContactInformation.Telem | * |
| DCM::ContentAuthorList | Projectgroep Generieke Overdrachtsgegevens & Kerngroep Registratie aan de Bron |
| DCM::CreationDate | 7-6-2012 |
| DCM::DeprecatedDate | |
| DCM::DescriptionLanguage | nl |
| DCM::EndorsingAuthority.Address | |
| DCM::EndorsingAuthority.Name | PM |
| DCM::EndorsingAuthority.Telem | |
| DCM::Id | 2.16.840.1.113883.2.4.3.11.60.40.3.13.1 |
| DCM::KeywordList | laboratorium uitslag, lab, laboratorium bepaling |
| DCM::LifecycleStatus | Final |
| DCM::ModelerList | Kerngroep Registratie aan de Bron |
| DCM::Name | nl.zorg.LaboratoriumUitslag |
| DCM::PublicationDate | 04-09-2017 |
| DCM::PublicationStatus | Prepublished |
| DCM::ReviewerList | Projectgroep Generieke Overdrachtsgegevens & Kerngroep Registratie aan de Bron |
| DCM::RevisionDate | 04-09-2017 |
| DCM::Superseeds | nl.zorg.OverdrachtLaboratoriumUitslag-v3.0 |
| DCM::Version | 4.0 |
| HCIM::PublicationLanguage | EN |

1.1 Revision History

Publicatieversie 1.0 (15-02-2013)

Publicatieversie 1.1 (01-07-2013)

Publicatieversie 1.2 (01-04-2015)

Bevat: ZIB-238, ZIB-239, ZIB-240, ZIB-241, ZIB-242, ZIB-243, ZIB-244, ZIB-245, ZIB-246, ZIB-353, ZIB-361, ZIB-367, ZIB-370.

Incl. algemene wijzigingsverzoeken:

ZIB-94, ZIB-154, ZIB-200, ZIB-201, ZIB-309, ZIB-324, ZIB-326.

Publicatieversie 1.2.1 (22-05-2015)

Bevat: ZIB-392.

Publicatieversie 1.2.2 (16-07-2015)

Bevat: ZIB-420.

Publicatieversie 3.0 (01-05-2016)

Bevat: ZIB-423, ZIB-453.

Publicatieversie 4.0 (04-09-2017)

Bevat: ZIB-479, ZIB-549, ZIB-564, ZIB-576, ZIB-481, ZIB-577.

1.2 Concept

A laboratory result describes the result of a laboratory analysis.

These are specimen-oriented tests as performed in laboratories such as Clinical Chemistry, Serology, Microbiology, etc.

In addition to the results of tests with a singular result, this concept can also contain the results of more complex tests with multiple results or a ‘panel’.

1.3 Mindmap

1.4 Purpose

Laboratory tests are done for the purpose of diagnosing and preventing disease and follow-up on the effects of treatment.

1.5 Patient Population

1.6 Evidence Base


There are two information models for recording laboratory test results: TextResult and LaboratoryTestResult.

In the case of laboratory test results, it is difficult to clearly indicate exactly when to use this information model and when to use the TextResult information model.

In general, laboratory tests resulting in a value (7.1 mmol/L), ordinal number (++ from series to +++) or a quantitative result (Low) are recorded using this information model. The TextResult information model is better suited for textual results that are more descriptive in nature and which are longer than just a few words. Both types of tests occur in almost all laboratories.

The applicability of the aforementioned information models is not determined by the kind of lab but by the kind of result.

1.7 Information Model



| «rootconcept» | LaboratoryTestResult | |
|----------------|---|--|
| Definitie | Root concept of the LaboratoryTestResult information model. This root concept contains all data elements of the LaboratoryTestResult information model. | |
| Datatype | | |
| DCM::ConceptId | NL-CM:13.1.1 | |
| Opties | | |

| «data» | PanelOrBattery | |
|-------------------|---|---|
| Definitie | For laboratory tests comprising multiple subtests and often requested together as a whole, this concept contains the name of the compound request (often indicated as a ‘panel’, ‘battery’ or ‘cluster’). Examples include: blood gases and EBV serology. | |
| Datatype | CD | |
| DCM::ConceptId | NL-CM:13.1.4 | |
| DCM::ExampleValue | Bloedgassen | |
| DCM::ValueSet | PanelOrBatteryCodelist | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.5 |
| Opties | | |

| «data» | ResultStatus | |
|-------------------|---|---|
| Definitie | The status of the laboratory test result .If the laboratory test is an panel/cluster, this status reflects the status of the whole panel/cluster. If the status item per subtest is used too, this status must be in accordance with these status values. | |
| Datatype | CD | |
| DCM::ConceptId | NL-CM:13.1.6 | |
| DCM::ExampleValue | Definitief | |
| DCM::ValueSet | ResultStatusCodelist | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.8 |

| | |
|----------------------------|--|
| Opties | |
| «data» | Comment |
| Definitie | Comments, such as a textual interpretation or advice accompanying the result, for example. |
| Datatype | ST |
| DCM::ConceptId | NL-CM:13.1.5 |
| DCM::DefinitionCode | LOINC: 48767-8 Annotation comment |
| Opties | |

| | | |
|--------------------------|--|---|
| «data» | ResultType | |
| Definitie | The type of result defines the laboratory specialty under which the test is categorized. | |
| Datatype | CD | |
| DCM::ConceptId | NL-CM:13.1.7 | |
| DCM::ExampleValue | Klinische Chemie | |
| DCM::ValueSet | ResultTypeCodelist | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.1 |
| Opties | | |

| | | |
|---------------------------------|---|---|
| «data» | RelatedResult::LaboratoryTestResult | |
| Definitie | Reference to related tests, e.g. paired tests or sequential tests like gram staining and microbiological cultures | |
| Datatype | | |
| DCM::ConceptId | NL-CM:13.1.33 | |
| DCM::ReferencedConceptId | NL-CM:13.1.1 | This is a reference to the rootconcept of information model LaboratoryTestResult. |
| Opties | | |

| | | |
|---------------------------|---|-----------------------|
| «document» | PanelOrBatteryCodelist | |
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.5 | |
| Opties | | |
| OnderzoekCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.5 | |
| Codes | Coding Syst. Name | Coding System OID |
| Alle waarden | LOINC | 2.16.840.1.113883.6.1 |


Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

| | | | | |
|---------------------------------|---|---------------|----------------|-------------|
| «document» | ResultStatusCodelist | | | |
| Definitie | | | | |
| Datatype | | | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.8 | | | |
| Opties | | | | |
| ResultaatStatusCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.8 | | | |
| Concept Name | Concept Code | CodeSys. Name | CodeSystem OID | Description |

| | | | | |
|-------------|-------------|-----------------|---|---------------|
| Pending | pending | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Uitslag volgt |
| Preliminary | preliminary | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Voorlopig |
| Final | final | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Definitief |
| Appended | appended | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Aanvullend |
| Corrected | corrected | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Gecorrigeerd |

| «document» | ResultTypeCodelist | | | |
|-------------------------------|---|-------------------|---|------------------------|
| Definitie | | | | |
| Datatype | | | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.1 | | | |
| Opties | | | | |
| ResultaatTypeCodelijst | | | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.1 | |
| Concept Name | Concept Code | Coding Syst. Name | Coding System OID | Omschrijving |
| Hematology | 252275004 | SNOMED CT | 2.16.840.1.113883.6.96 | Hematologie |
| Chemistry | 275711006 | SNOMED CT | 2.16.840.1.113883.6.96 | Klinische chemie |
| Serology | 68793005 | SNOMED CT | 2.16.840.1.113883.6.96 | Serologie/ immunologie |
| Virology | 395124008 | SNOMED CT | 2.16.840.1.113883.6.96 | Virologie |
| Toxicology | 314076009 | SNOMED CT | 2.16.840.1.113883.6.96 | Toxicologie |
| Microbiology | 19851009 | SNOMED CT | 2.16.840.1.113883.6.96 | Microbiologie |
| Molecular genetics | 405825005 | SNOMED CT | 2.16.840.1.113883.6.96 | Moleculaire genetica |

1.7.1 LaboratoriumTest



| «container» | LaboratoryTest |
|----------------|---|
| Definitie | Container of the LaboratoryTest concept. This container contains all data elements of the LaboratoryTest concept. |
| Datatype | |
| DCM::ConceptId | NL-CM:13.1.3 |
| Opties | |

| «data» | Test |
|-------------------|---|
| Definitie | The name and code of the executed test. |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.8 |
| DCM::ExampleValue | HbA1c |
| DCM::ValueSet | TestNameCodelist OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.3 |
| Opties | |

| «data» | TestMethod |
|---------------------|---|
| Definitie | The test method used to obtain the result. |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.9 |
| DCM::DefinitionCode | SNOMED CT: 246501002 Technique (attribute) |
| DCM::ExampleValue | IFCC |
| DCM::ValueSet | TestMethodCodelist OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.4 |

| | |
|---------------|--|
| Opties | |
|---------------|--|

| | |
|--------------------------|--|
| «data» | TestDateTime |
| Definitie | The date and if possible the time at which the test was carried out. |
| Datatype | TS |
| DCM::ConceptId | NL-CM:13.1.13 |
| DCM::ExampleValue | 10-07-2012 20:15 |
| Opties | |

| | |
|--------------------------|---|
| «data» | Result |
| Definitie | The test result. Depending on the type of test, the result will consist of a value with a unit or a coded value (ordinal or nominal). |
| Datatype | ANY |
| DCM::ConceptId | NL-CM:13.1.10 |
| DCM::ExampleValue | 53 mmol/mol |
| Opties | |

| | |
|-----------------------|---|
| «data» | TestResultStatus |
| Definitie | The status of the test result of this test. If the laboratory test is an panel/cluster, the overall status is given in the status of the panel/cluster. |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.31 |
| DCM::ValueSet | TestResultStatusCodelist |
| | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.15 |
| Opties | |

| | |
|--------------------------|--|
| «data» | ReferenceRangeUpperLimit |
| Definitie | The upper reference limit for the patient of the value measured in the test. |
| Datatype | ANY |
| DCM::ConceptId | NL-CM:13.1.11 |
| DCM::ExampleValue | 42 mmol/mol |
| Opties | |

| | |
|--------------------------|--|
| «data» | ReferenceRangeLowerLimit |
| Definitie | The lower reference limit for the patient of the value measured with the test. |
| Datatype | ANY |
| DCM::ConceptId | NL-CM:13.1.12 |
| DCM::ExampleValue | 20 mmol/mol |
| Opties | |

| | |
|-----------------------|---|
| «data» | InterpretationMethod |
| Definitie | The method used to determine interpretation flags. An example of this is EUCAST, for determining clinical breakpoints in microbiological susceptibility tests |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.30 |

| | | |
|---------------|------------------------------|--|
| DCM::ValueSet | InterpretationMethodCodelist | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.14 |
| Opties | | |

| «data» | ResultFlags | |
|---------------------|---|---|
| Definitie | Attention codes indicating whether the result is above or below certain reference values or interpreting the result otherwise.(Resistant) | |
| Datatype | CD | |
| DCM::ConceptId | NL-CM:13.1.14 | |
| DCM::DefinitionCode | SNOMED CT: 363713009 Has interpretation | |
| DCM::ExampleValue | High | |
| DCM::ValueSet | ResultFlagsCodelist | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.7 |
| Opties | | |

| «data» | ResultInterpretation | |
|---------------------|--|--|
| Definitie | Comment of the laboratory specialist regarding the interpretation of the results | |
| Datatype | ST | |
| DCM::ConceptId | NL-CM:13.1.32 | |
| DCM::DefinitionCode | SNOMED CT: 441742003 Evaluation finding | |
| Opties | | |

| «document» | TestNameCodelist | |
|-----------------|--|-----------------------|
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.3 | |
| Opties | | |
| TestCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.3 | |
| Codes | Coding Syst. Name | Coding System OID |
| Alle waarden | LOINC | 2.16.840.1.113883.6.1 |

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

| «document» | TestMethodCodelist | |
|--|--|------------------------|
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.4 | |
| Opties | | |
| TestmethodeCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.4 | |
| Codes | Coding Syst. Name | Coding System OID |
| SNOMED CT: < 272394005 Technique (qualifier value) | SNOMED CT | 2.16.840.1.113883.6.96 |

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

| «document» | TestResultStatusCodelist | |
|------------|--------------------------|--|
| Opties | | |

| Definitie | | | | |
|-----------------------------------|--|--|---|-----------------------------------|
| Datatype | | | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.15 | | | |
| Opties | | | | |
| TestUitslagStatusCodelijst | | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.15 | | |
| Concept Name | Concept Code | CodeSys. Name | CodeSystem OID | Description |
| Pending | pending | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Uitslag volgt .11.60.40.4.16.1 |
| Preliminary | preliminary | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Voorlopig .11.60.40.4.16.1 |
| Final | final | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Definitief .11.60.40.4.16.1 |
| Appended | appended | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Aanvullend .11.60.40.4.16.1 |
| Corrected | corrected | ResultaatStatus | 2.16.840.1.113883.2.4.3 .11.60.40.4.16.1 | Gecorrigeerd .11.60.40.4.16.1 |


| | | | | |
|------------------------|--|--|--|--|
| «document» | InterpretationMethodCodelist | | | |
| Definitie | | | | |
| Datatype | | | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.14 | | | |
| Opties | | | | |

| InterpretatieMethodeCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.14 | | | |
|--------------------------------------|--|-----------------|------------------------|-------------|
| Concept Name | Concept Code | CodeSystem Name | CodeSystem OID | Description |
| EUCAST | tbd | SNOMED CT | 2.16.840.1.113883.6.96 | EUCAST |

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

| «document» | ResultFlagsCodelist | | | |
|--------------------------------------|---|---|------------------------|------------------------|
| Definitie | | | | |
| Datatype | | | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.7 | | | |
| Opties | | | | |
| InterpretatieVlaggenCodelijst | | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.7 | | |
| Concept Name | Concept Code | Coding Syst. Name | Coding System OID | Description |
| Above reference range | 281302008 | SNOMED CT | 2.16.840.1.113883.6.96 | Boven referentiewaarde |
| Below reference range | 281300000 | SNOMED CT | 2.16.840.1.113883.6.96 | Onder referentiewaarde |
| Intermediate | 11896004 | SNOMED CT | 2.16.840.1.113883.6.96 | Intermediair |
| Resistant | 30714006 | SNOMED CT | 2.16.840.1.113883.6.96 | Resistent |
| Susceptible | 131196009 | SNOMED CT | 2.16.840.1.113883.6.96 | Sensitief |

1.7.2 Monster



| «container» | Specimen |
|----------------------------|---|
| Definitie | Container of the Specimen concept. This container contains all data elements of the Specimen concept. |
| Datatype | |
| DCM::ConceptId | NL-CM:13.1.2 |
| DCM::DefinitionCode | SNOMED CT: 123038009 monster |
| Opties | |

| «data» | SpecimenId |
|-----------------------|---|
| Definitie | Identification number of the material obtained, as a reference for inquiries to the source organization. In a transmural setting, this number will consist of a specimen number including the identification of the issuing organization, to be unique outside of the borders of an organization. |
| Datatype | II |
| DCM::ConceptId | NL-CM:13.1.15 |
| Opties | |

| «data» | SpecimenNumberExtension |
|------------------|--|
| Definitie | The specimen number extension is used when the collected material is distributed from the original tube or container across multiple tubes. In combination with the specimen Id the extension yield a unique |

| | |
|-----------------------|---|
| | identification of the tube or container |
| Datatype | INT |
| DCM::ConceptId | NL-CM:13.1.20 |
| Opties | |

| «data» | ContainerType |
|-----------------------|---|
| Definitie | Container type describes the envelope in which the material is collected or sent. Examples include blood tubes, transport container, possibly including culture medium. |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.21 |
| DCM::ValueSet | ContainerTypeCodelist |
| | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.9 |
| Opties | |

| «data» | SpecimenMaterial |
|----------------------------|--|
| Definitie | SpecimenMaterial describes the material obtained. If the LOINC test code also implicitly describes a material, this element may not conflict with the description. If desired, this component can provide a more detailed description of the material: LOINC codes only contain the materials at a main level. This is in line with the agreements made in the IHE/Nictiz program e-Lab. If the test is carried out on derived material (such as plasma), this element will still contain the material drawn (in this case, blood). In this case, the LOINC code will generally refer to plasma. |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.16 |
| DCM::DefinitionCode | SNOMED CT: 370133003 Specimen substance |
| DCM::ExampleValue | Urine |
| DCM::ValueSet | SpecimenMaterialCodelist |
| | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.6 |
| Opties | |

| «data» | Microorganism |
|-----------------------|---|
| Definitie | In particular in microbiological determinations the subject of the test is an isolate of certain microorganism rather than a material. This concept provides the ability to capture information about this microorganism. |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.22 |
| DCM::ValueSet | MicroorganismCodelist |
| | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.10 |
| Opties | |

| «data» | CollectedVolume |
|------------------|--|
| Definitie | Total volume of the collected material. If it is necessary to determine the absolute amount of a particular substance in the collected material, the volume thereof must be given. |
| Datatype | PQ |

| | | |
|-----------------------|---------------|--|
| DCM::ConceptId | NL-CM:13.1.23 | |
| Opties | | |

| «data» | CollectionPeriod |
|-----------------------|---|
| Definitie | If the material has not been collected at a single point in time but over a certain period, this period can be captured in this concept. An example is 24 hour urine. |
| Datatype | TS |
| DCM::ConceptId | NL-CM:13.1.24 |
| Opties | |

| «data» | CollectionDateTime |
|----------------------------|--|
| Definitie | Time at which the material was collected. |
| Datatype | TS |
| DCM::ConceptId | NL-CM:13.1.17 |
| DCM::DefinitionCode | SNOMED CT: 399445004 specimen collection date |
| DCM::ExampleValue | 10-07-2012 17:20 |
| Opties | |

| «data» | ??DateTime |
|-----------------------|---|
| Definitie | Date and time that the material is handed over at the laboratory or specimen collection center. This is the issue with material that is collected by the patient himself. |
| Datatype | TS |
| DCM::ConceptId | NL-CM:13.1.25 |
| Opties | |

| «data» | CollectionMethod |
|----------------------------|---|
| Definitie | If relevant for the results, the method of obtaining the specimen can be entered as well. |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.18 |
| DCM::DefinitionCode | SNOMED CT: 118171006 Specimen procedure |
| DCM::ExampleValue | Midstream |
| DCM::ValueSet | CollectioMethodCodelist |
| | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.2 |
| Opties | |

| «data» | AnatomicalLocation |
|----------------------------|---|
| Definitie | Anatomic location where the material is collected, e.g. elbow |
| Datatype | CD |
| DCM::ConceptId | NL-CM:13.1.26 |
| DCM::DefinitionCode | SNOMED CT: 405814001 Procedure site - Indirect |
| DCM::ValueSet | SpecimenAnatomicalLocatio nCodelist |
| | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.11 |
| Opties | |

| «data» | Laterality | |
|---------------------|---|--|
| Definitie | Laterality adds information about body side to the anatomic location, e.g. left | |
| Datatype | CD | |
| DCM::ConceptId | NL-CM:13.1.27 | |
| DCM::DefinitionCode | SNOMED CT: 272741003 Laterality | |
| DCM::ValueSet | LateralityCodelist | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.12 |
| Opties | | |

| «data» | Morphology | |
|---------------------|--|--|
| Definitie | Morphology describes morphological abnormalities of the anatomical location where the material is taken, for example wound, ulcer. | |
| Datatype | CD | |
| DCM::ConceptId | NL-CM:13.1.28 | |
| DCM::DefinitionCode | SNOMED CT: 118168003 Specimen source morphology | |
| DCM::ValueSet | MorphologyCodelist | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.13 |
| Opties | | |

| «data» | SpecimenSource | |
|---------------------|---|--|
| Definitie | If the material is not collected directly from the patient but comes from a patient-related object, e.g. a cathetertip, this source of material can be recorded here. | |
| Datatype | ST | |
| DCM::ConceptId | NL-CM:13.1.29 | |
| DCM::DefinitionCode | SNOMED CT: 127454002 Device specimen | |
| Opties | | |

| «data» | Comment | |
|---------------------|---|--|
| Definitie | Comments on administering the test, such as drawing material after a (glucose) stimulus or taking medicine. | |
| Datatype | ST | |
| DCM::ConceptId | NL-CM:13.1.19 | |
| DCM::DefinitionCode | LOINC: 48767-8 Annotation comment | |
| DCM::ExampleValue | Na (glucose)stimulus | |
| Opties | | |

| «document» | SpecimenMaterialCodelist | |
|-----------------|---|--|
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.6 | |
| Opties | | |

| Codes | Coding Syst. Name | Coding System OID |
|-----------------------------------|-------------------|------------------------|
| SNOMED CT: <105590001 substantie | SNOMED CT | 2.16.840.1.113883.6.96 |

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

| «document» | ContainerTypeCodelist | |
|---|--|------------------------|
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.9 | |
| Opties | | |
| ContainerTypeCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.9 | |
| Codes | Coding Syst. Name | Coding System OID |
| SNOMED CT: <434711009 Specimen container (physical object) | SNOMED CT | 2.16.840.1.113883.6.96 |

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

| «document» | MicroorganismCodelist | |
|---|---|------------------------|
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.10 | |
| Opties | | |
| MicroorganismeCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.10 | |
| Codes | Coding Syst. Name | Coding System OID |
| SNOMED CT: ^2581000146104 simple referentieset voor micro-organismen (foundation metadata concept) | SNOMED CT | 2.16.840.1.113883.6.96 |

| «document» | CollectioMethodCodelist | |
|---|--|------------------------|
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.2 | |
| Opties | | |
| AfnameprocedureCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.2 | |
| Codes | Coding Syst. Name | Coding System OID |
| SNOMED CT: <17636008 specimen collection | SNOMED CT | 2.16.840.1.113883.6.96 |

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

| «document» | SpecimenAnatomicalLocationCodelist | |
|------------------------------------|---|-------------------|
| Definitie | | |
| Datatype | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.11 | |
| Opties | | |
| MonsterAnatomischeLocatieCodelijst | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.11 | |
| Codes | Coding Syst. Name | Coding System OID |

| | | |
|---|-----------|------------------------|
| SNOMED CT: < 442083009 Anatomical or acquired body structure | SNOMED CT | 2.16.840.1.113883.6.96 |
|---|-----------|------------------------|

| «document» | LateralityCodelist | | | |
|-----------------------|--|---|------------------------|-----------------|
| Definitie | | | | |
| Datatype | | | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.12 | | | |
| Opties | | | | |
| LateraliteitCodelijst | | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.12 | | |
| Concept Name | Concept Code | CodeSys. Name | CodeSystem OID | Description |
| Left | 7771000 | SNOMED CT | 2.16.840.1.113883.6.96 | Links |
| Right | 24028007 | SNOMED CT | 2.16.840.1.113883.6.96 | Rechts |
| Right and left | 51440002 | SNOMED CT | 2.16.840.1.113883.6.96 | Rechts en links |

| «document» | MorphologyCodelist | | | |
|---|--|---|--|--|
| Definitie | | | | |
| Datatype | | | | |
| DCM::ValueSetId | 2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.13 | | | |
| Opties | | | | |
| MorfologieCodelijst | | OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.13 | | |
| Codes | Coding Syst. Name | Coding System OID | | |
| SNOMED CT: < 49755003 Morphologically abnormal structure | SNOMED CT | 2.16.840.1.113883.6.96 | | |

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

1.8 Example Instances

| LaboratoriumUitslag | | | | | | | | | | |
|---------------------|------------------|-------------------|------------------|------------------|------------------|------------|-----------------------|-----------------------|-------------------|--|
| Resultaat Type | Resultaat Status | Monster | | LaboratoriumTest | | | | | | |
| | | Monster materiaal | Afname DatumTijd | TestNaam | Test DatumTijd | Uitslag | Referentie Ondergrens | Referentie Bovengrens | Resultaat Vlaggen | |
| Klinische chemie | Definitief | Bloed | 12-06-2012 09:00 | Natrium | 12-06-2012 13:15 | 138 mmol/l | 136 mmol/l | 146 mmol/l | | |

| LaboratoriumUitslag | | | | | | | | | | |
|---------------------|------------------|-------------------|------------------|------------------|------------------|------------|-----------------------|-----------------------|------------------------|--|
| Resultaat Type | Resultaat Status | Monster | | LaboratoriumTest | | | | | | |
| | | Monster materiaal | Afname DatumTijd | TestNaam | Test DatumTijd | Uitslag | Referentie Ondergrens | Referentie Bovengrens | Resultaat Vlaggen | |
| Klinische chemie | Definitief | Bloed | 23-05-2012 08:08 | Chloride | 23-05-2012 12:00 | 109 mmol/l | 99 mmol/l | 108 mmol/l | Boven referentiewaarde | |

| LaboratoriumUitslag | | | | | | | | | | |
|---------------------|------------------|-------------------|------------------|------------------|------------------|----------|-----------------------|-----------------------|-------------------|--|
| Resultaat Type | Resultaat Status | Monster | | LaboratoriumTest | | | | | | |
| | | Monster materiaal | Afname DatumTijd | TestNaam | Test DatumTijd | Uitslag | Referentie Ondergrens | Referentie Bovengrens | Resultaat Vlaggen | |
| Virologie | Definitief | Bloed | 16-01-2012 08:00 | Hepatitis A IgM | 16-01-2012 10:12 | Negatief | | | | |

1.9 Instructions

1.10 Interpretation

1.11 Care Process

1.12 Example of the Instrument

1.13 Constraints

1.14 Issues

1.15 References

1. Nederlandse Vereniging voor Medische Microbiologie (2010) *ELab en EvT*. [Online] Beschikbaar op: http://www.nvmm.nl/ict/vereniging/werkgroepen_commissies/elab-en-evt [Geraadpleegd: 23 juli 2014].

1.16 Functional Model

1.17 Traceability to other Standards

1.18 Disclaimer

This Health and Care Information Model (a.k.a Clinical Building Block) has been made in collaboration with several different parties in healthcare. These parties asked Nictiz to manage good maintenance and development of the information models. Hereafter, these parties and Nictiz are referred to as the collaborating parties. The collaborating parties paid utmost attention to the reliability and topicality of the data in these Health and Care Information Models. Omissions and inaccuracies may however occur. The collaborating parties are not liable for any damages resulting from omissions or inaccuracies in the information provided, nor are they liable for damages resulting from problems caused by or inherent to distributing information on the internet, such as malfunctions, interruptions, errors or delays in information or services provide by the parties to you or by you to the parties via a website or via e-mail, or any other digital means. The collaborating parties will also not accept liability for any damages resulting from the use of data, advice or ideas provided by or on behalf of the parties by means of this Health and Care Information Model. The parties will not accept any liability for the content of information in this Health and Care Information Model to which or from which a hyperlink is referred. In the event of contradictions in mentioned Health and Care Information Model documents and files, the most recent and highest version of the listed order in the revisions will indicate the priority of the documents in question. If information included in the digital version of this Health and Care Information Model is also distributed in writing, the written version will be leading in case of textual differences. This will apply if both have the same version number and date. A definitive version has priority over a draft version. A revised version has priority over previous versions.

1.19 Terms of Use

The user may use the information in this Health and Care Information Model without limitations. The copyright provisions in the paragraph concerned apply to copying, distributing and passing on information from this Health and Care Information Model.

1.20 Copyrights

The user may copy, distribute and pass on the information in this Health and Care Information Model under the conditions that apply for Creative Commons license Attribution-NonCommercial-ShareAlike 3.0 Netherlands (CC BY-NCSA-3.0). The content is available under Creative Commons Attribution-NonCommercial-ShareAlike 3.0 (see also <http://creativecommons.org/licenses/by-nc-sa/3.0/nl/>)