

# Health&Care Information Model: nl.zorg.Blaasfunctie

Final

Managed by:

Better healthcare  
through better IT



# Content

<b>1. nl.zorg.Blaasfunctie-v3.0</b> .....	<b>3</b>
1.1 Revision History .....	3
1.2 Concept .....	3
1.3 Mindmap .....	3
1.4 Purpose .....	3
1.5 Patient Population .....	4
1.6 Evidence Base .....	4
1.7 Information Model .....	4
1.8 Example Instances .....	7
1.9 Instructions .....	8
1.10 Interpretation .....	8
1.11 Care Process .....	8
1.12 Example of the Instrument .....	8
1.13 Constraints .....	8
1.14 Issues .....	9
1.15 References .....	9
1.16 Functional Model .....	9
1.17 Traceability to other Standards .....	9
1.18 Disclaimer .....	9
1.19 Terms of Use .....	9
1.20 Copyrights .....	9

# 1. nl.zorg.Blaasfunctie-v3.0

DCM::CoderList	Werkgroep RadB Verpleegkundige Gegevens
DCM::ContactInformation.Address	*
DCM::ContactInformation.Name	*
DCM::ContactInformation.Telecom	*
DCM::ContentAuthorList	Werkgroep RadB Verpleegkundige Gegevens
DCM::CreationDate	1-4-2014
DCM::DeprecatedDate	
DCM::DescriptionLanguage	nl
DCM::EndorsingAuthority.Address	
DCM::EndorsingAuthority.Name	PM
DCM::EndorsingAuthority.Telecom	
DCM::Id	2.16.840.1.113883.2.4.3.11.60.40.3.4.14
DCM::KeywordList	Mictie, continentie
DCM::LifecycleStatus	Final
DCM::ModelerList	Werkgroep RadB Verpleegkundige Gegevens
DCM::Name	nl.zorg.Blaasfunctie
DCM::PublicationDate	1-5-2016
DCM::PublicationStatus	Published
DCM::ReviewerList	Projectgroep RadB Verpleegkundige Gegevens & Kerngroep Registratie aan de Bron
DCM::RevisionDate	8-9-2015
DCM::Superseeds	nl.nfu.Blaasfunctie-v1.0
DCM::Version	3.0

## 1.1 Revision History

Publicatieversie 1.0 (01-07-2015)

Publicatieversie 3.0 (01-05-2016)

Bevat: ZIB-453

## 1.2 Concept

**! Attention: this information model is undergoing major revision. A new version will be available after summer.**

The primary function of the bladder is temporary storage of urine. A second important function of the bladder is secretion of the stored urine at the moment there is cause to do so. These functions can be disrupted due to various causes.

## 1.3 Mindmap

## 1.4 Purpose

Information on bladder function disorders is important in drawing up the nursing plan and for patient care.


In a transfer situation, it offers the ability to anticipate the arrival of the patient and to achieve continuity in healthcare for the patient. This way, material or expertise can be coordinated beforehand, for example.

## 1.5 Patient Population

## 1.6 Evidence Base

The definitions of the concepts were (partially) taken from the ICF definitions.

## 1.7 Information Model



«rootconcept»	Blaasfunctie	
Alias	EN: BladderFunction	
Definition	Root concept of the BladderFunction building block. This root concept contains all data elements of the BladderFunction building block.	
Datatype		
DCM::DefinitionCode	SNOMED CT: 28278009	Micturition
DCM::DefinitionCode	NL-CM:4.14.1	
Options		

«data»	UrineContinentie	
Alias	EN: UrineContinece	
Definition	Urine continence describes control over urine discharge. This includes	

	disorders such as stress incontinence, urge incontinence, reflex incontinence, constant incontinence and mixed incontinence.	
<b>Datatype</b>	CO	
<b>DCM::DefinitionCode</b>	NL-CM:4.14.2	
<b>DCM::DefinitionCode</b>	ICF: b6202 Urinecontinentie	
<b>DCM::ExampleValue</b>	Lichte stoornis	
<b>DCM::ValueSet</b>	UrineContinentieCodelijst	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.4.14.1
<b>Options</b>		

<b>«data»</b>	<b>SoortIncontinentie::Probleem</b>	
<b>Alias</b>	EN: IncontinenceType::Problem	
<b>Definition</b>	A specification of the type of incontinence.	
<b>Datatype</b>		
<b>DCM::DefinitionCode</b>	SNOMED CT: 165232002 Urinary incontinence	
<b>DCM::DefinitionCode</b>	NL-CM:4.14.8	
<b>DCM::ExampleValue</b>	Stress incontinentie	
<b>DCM::ReferencedDefinitionCode</b>	NL-CM:5.1.2	This is a reference to concept Probleem in information model OverdrachtConcern.
<b>Options</b>		

<b>«data»</b>	<b>Urostoma::Stoma</b>	
<b>Alias</b>	EN: UrinaryStoma::Stoma	
<b>Definition</b>	UrinaryStoma describes the presence of an unnatural, surgically created opening. If the patient has a urinary stoma, the date of placement, location of the stoma and any stoma pouches and an explanation can be entered. The type of stoma is automatically a urinary stoma in this case.n	
<b>Datatype</b>		
<b>DCM::DefinitionCode</b>	NL-CM:4.14.3	
<b>DCM::DefinitionCode</b>	SNOMED CT: 281601009 Urological stoma	
<b>DCM::ReferencedDefinitionCode</b>	NL-CM:5.2.1	This is a reference to concept Stoma in information model Stoma.
<b>Options</b>		

<b>«data»</b>	<b>UrineKatheter::MedischHulpmiddel</b>	
<b>Alias</b>	EN: UrineCatheter::MedicalAid	
<b>Definition</b>	UrineCatheter describes the presence of a urine catheter. If this is the case, the date of placement, urine catheter location, explanation and catheter identification information can be described in addition to the type of catheter.	
<b>Datatype</b>		
<b>DCM::DefinitionCode</b>	SNOMED CT: 20568009 Urinary catheter	
<b>DCM::DefinitionCode</b>	NL-CM:4.14.4	
<b>DCM::ReferencedDefinitionCode</b>	NL-CM:10.1.1	This is a reference to concept MedischHulpmiddel in information model MedischHulpmiddel.
<b>Options</b>		

«data»	ProductType	
Alias	EN: ProductType	
Definition	A description of the type of catheter.	
Datatype	CD	
DCM::DefinitionCode	NL-CM:10.1.3	
DCM::ExampleValue	Verblijfskatheter	
DCM::ValueSet	UrineKatheterTypeCodelijst	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.4.14.2
Options		

«data»	Incontinentiemateriaal::MedischHulpmiddel	
Alias	EN: IncontinenceMaterial::MedicalAid	
Definition	A description of the incontinence material used.	
Datatype		
DCM::DefinitionCode	NL-CM:4.14.5	
DCM::DefinitionCode	SNOMED CT: 337588003 Incontinence appliance	
DCM::ExampleValue	Inlegverband	
DCM::ReferencedDefinitionCode	NL-CM:10.1.1	This is a reference to concept MedischHulpmiddel in information model MedischHulpmiddel.
Options		

«data»	Toelichting	
Alias	EN: Explanation	
Definition	Explanation of the urinary (in)continence and treatment.	
Datatype	ST	
DCM::DefinitionCode	NL-CM:4.14.6	
DCM::DefinitionCode	LOINC: 48767-8 Annotation comment	
Options		

«data»	VerpleegkundigeActie	
Alias	EN: NursingProcedure	
Definition	Description of the help needed for micturition and the required facilitation of any urine catheters.	
Datatype		
DCM::DefinitionCode	NL-CM:4.14.7	
DCM::ExampleValue	Ondersteuning bij verwisselen zakje stoma	
DCM::ReferencedDefinitionCode	NL-CM:14.2.9	This is a reference to concept VerpleegkundigeActie in information model VerpleegkundigeInterventie.
Options		

«data»	GeplandeVPKActie::GeplandeZorgActiviteit	
Alias	EN: PlannedNursingProcedure::PlannedCareActivity	
Definition	Planned (nursing) procedure or 'clinical reminder' for example in facilitating the catheter, in which information can be included such as the date on which the catheter must be replaced.	

<b>Datatype</b>		
<b>DCM::DefinitionCode</b>	NL-CM:4.14.9	
<b>DCM::ReferencedDefinitionCode</b>	NL-CM:16.1.1	This is a reference to concept GeplandeZorgActiviteit in information model OverdrachtGeplandeZorgActiviteit.
<b>Options</b>		

<b>«document»</b>	<b>UrineContinentieCodelijst</b>	
<b>Alias</b>		
<b>Definition</b>		
<b>Datatype</b>		
<b>Options</b>		

<b>UrineContinentieCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.4.14.1</b>		
Concept Name	Concept Code	CodeSys. Name	CodeSystem OID	Description
Geen stoornis	b6202.0	ICF	2.16.840.1.113883.6.254	Stoornis 0-4%
Lichte stoornis	b6202.1	ICF	2.16.840.1.113883.6.254	Stoornis 5-24%
Matige stoornis	b6202.2	ICF	2.16.840.1.113883.6.254	Stoornis 25-49%
Ernstige stoornis	b6202.3	ICF	2.16.840.1.113883.6.254	Stoornis 50-95%
Volledige stoornis	b6202.4	ICF	2.16.840.1.113883.6.254	Stoornis 96-100%

<b>«document»</b>	<b>UrineKatheterTypeCodelijst</b>	
<b>Alias</b>		
<b>Definition</b>		
<b>Datatype</b>		
<b>Options</b>		

<b>UrineKatheterTypeCodelijst</b>			<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.4.14.2</b>	
Concept Name	Concept Code	Coding Syst. Name	Coding System OID	Description
Indwelling urinary catheter	23973005	SNOMED CT	2.16.840.1.113883.6.96	Verblijfskatheter
Intermittent urethral drainage catheter, sterile	469715007	SNOMED CT	2.16.840.1.113883.6.96	Eenmalige katheter
Incontinence sheath	337636000	SNOMED CT	2.16.840.1.113883.6.96	Condoomkatheter
Suprapubic catheter	286861005	SNOMED CT	2.16.840.1.113883.6.96	Suprapubische katheter
Other	OTH	NullFlavor	2.16.840.1.113883.5.1008	Anders

## 1.8 Example Instances

Blaasfunctie (urostoma)	
UrineContinentie	-
SoortIncontinentie	-
Incontinentiemateriaal	-
Urostoma	
StomaType	Urostoma
AanlegDatum	12-01-2014
Toelichting	Extra aandacht voor peristomale huid
UrineKatheter	
ProductType	-
BeginDatum	-
Toelichting	-
VerpleegkundigeActie	
Activiteit	Ondersteunen bij verwisselen plak bij stoma.
Toelichting	Kan bij problemen contact opnemen met stomaverpleegkundige

Blaasfunctie (continentie stoornis)	
UrineContinentie	Matige stoornis
SoortIncontinentie	overloopcontinentie
Incontinentiemateriaal	-
Urostoma	
Sinds	-
Opmerking bij urostoma	-
UrineKatheter	
ProductType	Verblijfskatheter
BeginDatum	12-05-2014
Toelichting	-
VerpleegkundigeActie	
Activiteit	Ondersteunen bij het verwisselen van nachtopvangzak.
Toelichting	Partner verzorgt insteekopening van katheter

## 1.9 Instructions

## 1.10 Interpretation

## 1.11 Care Process

## 1.12 Example of the Instrument

## 1.13 Constraints



## 1.14 Issues

## 1.15 References

1. International Classification of Functioning Disability and Health (ICF) [Online] Beschikbaar op: <http://www.rivm.nl/who-fic/icf.htm> [Geraadpleegd: 13 februari 2015]

## 1.16 Functional Model

## 1.17 Traceability to other Standards

## 1.18 Disclaimer

This Clinical Building Block has been made in collaboration with several different parties in healthcare. These parties asked Nictiz to manage good maintenance and development of the information models. Hereafter, these parties and Nictiz are referred to as the collaborating parties. The collaborating parties paid utmost attention to the reliability and topicality of the data in these Health and Care Information Model (a.k.a Clinical Building Block)s. Omissions and inaccuracies may however occur. The collaborating parties are not liable for any damages resulting from omissions or inaccuracies in the information provided, nor are they liable for damages resulting from problems caused by or inherent to distributing information on the internet, such as malfunctions, interruptions, errors or delays in information or services provide by the parties to you or by you to the parties via a website or via e-mail, or any other digital means. The collaborating parties will also not accept liability for any damages resulting from the use of data, advice or ideas provided by or on behalf of the parties by means of this Clinical Building Block. The parties will not accept any liability for the content of information in this Clinical Building Block to which or from which a hyperlink is referred. In the event of contradictions in mentioned Clinical Building Block documents and files, the most recent and highest version of the listed order in the revisions will indicate the priority of the documents in question. If information included in the digital version of this Clinical Building Block is also distributed in writing, the written version will be leading in case of textual differences. This will apply if both have the same version number and date. A definitive version has priority over a draft version. A revised version has priority over previous versions.

## 1.19 Terms of Use

The user may use the information in this Clinical Building Block without limitations. The copyright provisions in the paragraph concerned apply to copying, distributing and passing on information from this Clinical Building Block.

## 1.20 Copyrights

The user may copy, distribute and pass on the information in this Clinical Building Block under the conditions that apply for Creative Commons license Attribution-NonCommercial-ShareAlike 3.0 Netherlands (CC BY-NC-SA-3.0). The content is available under Creative Commons Attribution-NonCommercial-ShareAlike 3.0 (see also <http://creativecommons.org/licenses/by-nc-sa/3.0/nl/>)