

# **Health & Care Information Model:**

## **nl.zorg.LaboratoryTestResult-v4.2**

Status:Final  
Release:2018  
Release status: Prepublished

Managed by:



# Content

1. nl.zorg.LaboratoryTestResult-v4.2 .....	3
1.1 Revision History .....	3
1.2 Concept .....	4
1.3 Mindmap .....	4
1.4 Purpose .....	4
1.5 Patient Population .....	4
1.6 Evidence Base .....	4
1.7 Information Model .....	4
1.7.1 LaboratoryTest .....	7
1.7.2 Specimen .....	12
1.8 Example Instances .....	17
1.9 Instructions .....	18
1.10 Interpretation .....	18
1.11 Care Process .....	18
1.12 Example of the Instrument .....	18
1.13 Constraints .....	18
1.14 Issues .....	18
1.15 References .....	18
1.16 Functional Model .....	18
1.17 Traceability to other Standards .....	18
1.18 Disclaimer .....	18
1.19 Terms of Use .....	19
1.20 Copyrights .....	19

# 1. nl.zorg.LaboratoryTestResult-v4.2

DCM::CoderList	Kerngroep Registratie aan de Bron
DCM::ContactInformation.Address	*
DCM::ContactInformation.Name	*
DCM::ContactInformation.Telco	*
DCM::ContentAuthorList	Projectgroep Generieke Overdrachtsgegevens & Kerngroep Registratie aan de Bron
DCM::CreationDate	7-6-2012
DCM::DeprecatedDate	
DCM::DescriptionLanguage	nl
DCM::EndorsingAuthority.Address	
DCM::EndorsingAuthority.Name	PM
DCM::EndorsingAuthority.Telco	
DCM::Id	2.16.840.1.113883.2.4.3.11.60.40.3.13.1
DCM::KeywordList	laboratorium uitslag, lab, laboratorium bepaling
DCM::LifecycleStatus	Final
DCM::ModelerList	Kerngroep Registratie aan de Bron
DCM::Name	nl.zorg.LaboratoriumUitslag
DCM::PublicationDate	01-10-2018
DCM::PublicationStatus	Prepublished
DCM::ReviewerList	Projectgroep Generieke Overdrachtsgegevens & Kerngroep Registratie aan de Bron
DCM::RevisionDate	31-12-2017
DCM::Superseeds	nl.zorg.LaboratoriumUitslag-v4.1
DCM::Version	4.2
HCIM::PublicationLanguage	EN

## 1.1 Revision History

Publicatieversie 1.0 (15-02-2013)

Publicatieversie 1.1 (01-07-2013)

Publicatieversie 1.2 (01-04-2015)

Bevat: ZIB-238, ZIB-239, ZIB-240, ZIB-241, ZIB-242, ZIB-243, ZIB-244, ZIB-245, ZIB-246, ZIB-353, ZIB-361, ZIB-367, ZIB-370.

Incl. algemene wijzigingsverzoeken:

ZIB-94, ZIB-154, ZIB-200, ZIB-201, ZIB-309, ZIB-324, ZIB-326.

Publicatieversie 1.2.1 (22-05-2015)

Bevat: ZIB-392.

Publicatieversie 1.2.2 (16-07-2015)

Bevat: ZIB-420.

Publicatieversie 3.0 (01-05-2016)

Bevat: ZIB-423, ZIB-453.

Publicatieversie 4.0 (04-09-2017)

Bevat: ZIB-479, ZIB-549, ZIB-564, ZIB-576, ZIB-481, ZIB-577.

Publicatieversie 4.1 (31-12-2017)

Bevat: ZIB-609, ZIB-611, ZIB-621, ZIB-645, ZIB-646.

## 1.2 Concept

A laboratory result describes the result of a laboratory analysis.

These are specimen-oriented tests as performed in laboratories such as Clinical Chemistry, Serology, Microbiology, etc.

In addition to the results of tests with a singular result, this concept can also contain the results of more complex tests with multiple results or a 'panel'.

## 1.3 Mindmap

## 1.4 Purpose

Laboratory tests are done for the purpose of diagnosing and preventing disease and follow-up on the effects of treatment.

## 1.5 Patient Population

## 1.6 Evidence Base


There are two information models for recording laboratory test results: TextResult and LaboratoryTestResult.

In the case of laboratory test results, it is difficult to clearly indicate exactly when to use this information model and when to use the TextResult information model.

In general, laboratory tests resulting in a value (7.1 mmol/L), ordinal number (++ from series to +++) or a quantitative result (Low) are recorded using this information model. The TextResult information model is better suited for textual results that are more descriptive in nature and which are longer than just a few words. Both types of tests occur in almost all laboratories.

The applicability of the aforementioned information models is not determined by the kind of lab but by the kind of result.

## 1.7 Information Model



<b>«rootconcept»</b> <b>LaboratoryTestResult</b>	
<b>Definitie</b>	Root concept of the LaboratoryTestResult information model. This root concept contains all data elements of the LaboratoryTestResult information model.
<b>Datatype</b>	
<b>DCM::ConceptId</b>	NL-CM:13.1.1
<b>Opties</b>	

<b>«data»</b> <b>PanelOrBattery</b>		
<b>Definitie</b>	For laboratory tests comprising multiple subtests and often requested together as a whole, this concept contains the name of the compound request (often indicated as a 'panel', 'battery' or 'cluster'). Examples include: blood gases and EBV serology.	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.4	
<b>DCM::ExampleValue</b>	Bloedgassen	
<b>DCM::ValueSet</b>	PanelOrBatteryCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.5
<b>Opties</b>		

<b>«data»</b> <b>ResultStatus</b>		
<b>Definitie</b>	The status of the laboratory test result .If the laboratory test is an panel/cluster, this status reflects the status of the whole panel/cluster. If the status item per subtest is used too, this status must be in accordance with these status values.	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.6	
<b>DCM::ExampleValue</b>	Definitief	
<b>DCM::ValueSet</b>	ResultStatusCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.8
<b>Opties</b>		

«data»	Comment
Definitie	Comments, such as a textual interpretation or advice accompanying the result, for example.
Datatype	ST
DCM::ConceptId	NL-CM:13.1.5
DCM::DefinitionCode	LOINC: 48767-8 Annotation comment
Opties	

«data»	ResultType	
Definitie	The type of result defines the laboratory specialty under which the test is categorized.	
Datatype	CD	
DCM::ConceptId	NL-CM:13.1.7	
DCM::ExampleValue	Klinische Chemie	
DCM::ValueSet	ResultTypeCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.1
Opties		

«data»	RelatedResult::LaboratoryTestResult	
Definitie	Reference to related tests, e.g. paired tests or sequential tests like gram staining and microbiological cultures	
Datatype		
DCM::ConceptId	NL-CM:13.1.33	
DCM::ReferencedConceptId	NL-CM:13.1.1	This is a reference to the rootconcept of information model LaboratoryTestResult.
Opties		

«context»	Requester::HealthProfessional	
Definitie	The healthcare provider and/or organization where or by whom the LaboratoryTestResult was requested.	
Datatype		
DCM::ConceptId	NL-CM:13.1.34	
DCM::ReferencedConceptId	NL-CM:17.1.1	This is a reference to the rootconcept of information model HealthProfessional.
Opties		

«context»	Performer::HealthcareProvider	
Definitie	The healthcare provider and/or organization where or by whom the LaboratoryTestResult was performed.	
Datatype		
DCM::ConceptId	NL-CM:13.1.35	
DCM::ReferencedConceptId	NL-CM:17.2.1	This is a reference to the rootconcept of information model HealthcareProvider.
Opties		

«document»	PanelOrBatteryCodelist
Definitie	


<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.5	
<b>Opties</b>		
<b>OnderzoekCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.5</b>
Codes	Coding Syst. Name	Coding System OID
Alle waarden	LOINC	2.16.840.1.113883.6.1

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

<b>«document»</b>	ResultStatusCodelist			
<b>Definitie</b>				
<b>Datatype</b>				
<b>DCM::ValueSetBinding</b>	Extensible			
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.8			
<b>Opties</b>				
<b>ResultaatStatusCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.8</b>		
Concept Name	Concept Code	CodeSys. Name	CodeSystem OID	Description
Pending	pending	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Uitslag volgt .11.60.40.4.16.1
Preliminary	preliminary	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Voorlopig .11.60.40.4.16.1
Final	final	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Definitief .11.60.40.4.16.1
Appended	appended	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Aanvullend .11.60.40.4.16.1
Corrected	corrected	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Gecorrigeerd .11.60.40.4.16.1

<b>«document»</b>	ResultTypeCodelist			
<b>Definitie</b>				
<b>Datatype</b>				
<b>DCM::ValueSetBinding</b>	Extensible			
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.1			
<b>Opties</b>				
<b>ResultaatTypeCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.1</b>		
Concept Name	Concept Code	Coding Syst. Name	Coding System OID	Omschrijving
Hematology	252275004	SNOMED CT	2.16.840.1.113883.6.96	Hematologie
Chemistry	275711006	SNOMED CT	2.16.840.1.113883.6.96	Klinische chemie
Serology	68793005	SNOMED CT	2.16.840.1.113883.6.96	Serologie/ immunologie
Virology	395124008	SNOMED CT	2.16.840.1.113883.6.96	Virologie
Toxicology	314076009	SNOMED CT	2.16.840.1.113883.6.96	Toxicologie
Microbiology	19851009	SNOMED CT	2.16.840.1.113883.6.96	Microbiologie
Molecular genetics	405825005	SNOMED CT	2.16.840.1.113883.6.96	Moleculaire genetica

## 1.7.1 LaboratoryTest



<b>«container»</b>	LaboratoryTest
<b>Definitie</b>	Container of the LaboratoryTest concept. This container contains all data elements of the LaboratoryTest concept.
<b>Datatype</b>	
<b>DCM::ConceptId</b>	NL-CM:13.1.3
<b>Opties</b>	

<b>«data»</b>	TestCode
<b>Definitie</b>	The name and code of the executed test.
<b>Datatype</b>	CD
<b>DCM::ConceptId</b>	NL-CM:13.1.8
<b>DCM::ExampleValue</b>	HbA1c
<b>DCM::ValueSet</b>	TestNameCodelist OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.3
<b>Opties</b>	

<b>«data»</b>	TestMethod
<b>Definitie</b>	The test method used to obtain the result.
<b>Datatype</b>	CD
<b>DCM::ConceptId</b>	NL-CM:13.1.9
<b>DCM::DefinitionCode</b>	SNOMED CT: 246501002 Technique (attribute)
<b>DCM::ExampleValue</b>	IFCC

<b>DCM::ValueSet</b>	TestMethodCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.4
<b>Opties</b>		

<b>«data»</b>	TestDateTime
<b>Definitie</b>	The date and if possible the time at which the test was carried out.
<b>Datatype</b>	TS
<b>DCM::ConceptId</b>	NL-CM:13.1.13
<b>DCM::ExampleValue</b>	10-07-2012 20:15
<b>Opties</b>	

<b>«data»</b>	TestResult
<b>Definitie</b>	The test result. Depending on the type of test, the result will consist of a value with a unit or a coded value (ordinal or nominal).
<b>Datatype</b>	ANY
<b>DCM::ConceptId</b>	NL-CM:13.1.10
<b>DCM::ExampleValue</b>	53 mmol/mol
<b>Opties</b>	

<b>«data»</b>	TestResultStatus
<b>Definitie</b>	The status of the test result of this test. If the laboratory test is an panel/cluster, the overall status is given in the status of the panel/cluster.
<b>Datatype</b>	CD
<b>DCM::ConceptId</b>	NL-CM:13.1.31
<b>DCM::ValueSet</b>	TestResultStatusCodelist
	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.15
<b>Opties</b>	

<b>«data»</b>	ReferenceRangeUpperLimit
<b>Definitie</b>	The upper reference limit for the patient of the value measured in the test.
<b>Datatype</b>	ANY
<b>DCM::ConceptId</b>	NL-CM:13.1.11
<b>DCM::ExampleValue</b>	42 mmol/mol
<b>Opties</b>	

<b>«data»</b>	ReferenceRangeLowerLimit
<b>Definitie</b>	The lower reference limit for the patient of the value measured with the test.
<b>Datatype</b>	ANY
<b>DCM::ConceptId</b>	NL-CM:13.1.12
<b>DCM::ExampleValue</b>	20 mmol/mol
<b>Opties</b>	

<b>«data»</b>	InterpretationMethod
<b>Definitie</b>	The method used to determine interpretation flags. An example of this is EUCAST, for determining clinical breakpoints in microbiological susceptibility tests
<b>Datatype</b>	CD

<b>DCM::ConceptId</b>	NL-CM:13.1.30	
<b>DCM::ValueSet</b>	InterpretationMethodCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.14
<b>Opties</b>		

«data»	ResultFlags	
<b>Definitie</b>	Attention codes indicating whether the result is above or below certain reference values or interpreting the result otherwise.(Resistant)	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.14	
<b>DCM::DefinitionCode</b>	SNOMED CT: 363713009 Has interpretation	
<b>DCM::ExampleValue</b>	High	
<b>DCM::ValueSet</b>	ResultFlagsCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.7
<b>Opties</b>		

«data»	ResultInterpretation	
<b>Definitie</b>	Comment of the laboratory specialist regarding the interpretation of the results	
<b>Datatype</b>	ST	
<b>DCM::ConceptId</b>	NL-CM:13.1.32	
<b>DCM::DefinitionCode</b>	SNOMED CT: 441742003 Evaluation finding	
<b>Opties</b>		

«document»	TestNameCodelist	
<b>Definitie</b>		
<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11.60.40.2.13.1.3	
<b>Opties</b>		
<b>TestCodeCodelijst</b>	<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.3</b>	
Codes	Coding Syst. Name	Coding System OID
Alle waarden	LOINC	2.16.840.1.113883.6.1

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

«document»	TestMethodCodelist	
<b>Definitie</b>		
<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11.60.40.2.13.1.4	
<b>Opties</b>		
<b>TestmethodeCodelijst</b>	<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.4</b>	
Codes	Coding Syst. Name	Coding System OID
SNOMED CT: < 272394005	SNOMED CT	2.16.840.1.113883.6.96  Technique (qualifier value)

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

<b>«document»</b>	TestResultStatusCodelist			
<b>Definitie</b>				
<b>Datatype</b>				
<b>DCM::ValueSetBinding</b>	Extensible			
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.15			
<b>Opties</b>				
<b>TestUitslagStatusCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.15</b>		
Concept Name	Concept Code	CodeSys. Name	CodeSystem OID	Description
Pending	pending	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Uitslag volgt .11.60.40.4.16.1
Preliminary	preliminary	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Voorlopig .11.60.40.4.16.1
Final	final	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Definitief .11.60.40.4.16.1
Appended	appended	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Aanvullend .11.60.40.4.16.1
Corrected	corrected	ResultaatStatus	2.16.840.1.113883.2.4.3 .11.60.40.4.16.1	Gecorrigeerd .11.60.40.4.16.1


<b>«document»</b>	InterpretationMethodCodelist			
<b>Definitie</b>				
<b>Datatype</b>				
<b>DCM::ValueSetBinding</b>	Extensible			
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.14			
<b>Opties</b>				
<b>InterpretatieMethodeCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.14</b>		
Concept Name	Concept Code	CodeSystem Name	CodeSystem OID	Description
EUCAST	tbd	SNOMED CT	2.16.840.1.113883.6.96	EUCAST

Specifiekeren waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

<b>«document»</b>	ResultFlagsCodelist			
<b>Definitie</b>				
<b>Datatype</b>				
<b>DCM::ValueSetBinding</b>	Extensible			
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.7			
<b>Opties</b>				
<b>InterpretatieVlaggenCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.7</b>		
Concept Name	Concept Code	Coding Syst. Name	Coding System OID	Description
Above reference range	281302008	SNOMED CT	2.16.840.1.113883.6.96	Boven referentiewaarde
Below reference range	281300000	SNOMED CT	2.16.840.1.113883.6.96	Onder referentiewaarde

Intermediate	11896004	SNOMED CT	2.16.840.1.113883.6.96	Intermediair
Resistant	30714006	SNOMED CT	2.16.840.1.113883.6.96	Resistent
Susceptible	131196009	SNOMED CT	2.16.840.1.113883.6.96	Sensitief

## 1.7.2 Specimen



«container»	Specimen
<b>Definitie</b>	Container of the Specimen concept. This container contains all data elements of the Specimen concept.
<b>Datatype</b>	
<b>DCM::ConceptId</b>	NL-CM:13.1.2
<b>DCM::DefinitionCode</b>	SNOMED CT: 123038009 monster
<b>Opties</b>	

«data»	SpecimenId
<b>Definitie</b>	Identification number of the material obtained, as a reference for inquiries to the source organization. In a transmural setting, this number will consist of a specimen number including the identification of the issuing organization, to be unique outside of the borders of an organization.
<b>Datatype</b>	II
<b>DCM::ConceptId</b>	NL-CM:13.1.15

<b>Opties</b>	
---------------	--

<b>«data»</b>	SpecimenNumberExtension
<b>Definitie</b>	The specimen number extension is used when the collected material is distributed from the original tube or container across multiple tubes. In combination with the specimen Id the extension yield a unique identification of the tube or container
<b>Datatype</b>	INT
<b>DCM::ConceptId</b>	NL-CM:13.1.20
<b>Opties</b>	

<b>«data»</b>	ContainerType	
<b>Definitie</b>	Container type describes the envelope in which the material is collected or sent. Examples include blood tubes, transport container, possibly including culture medium.	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.21	
<b>DCM::ValueSet</b>	ContainerTypeCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.9
<b>Opties</b>		

<b>«data»</b>	SpecimenMaterial	
<b>Definitie</b>	SpecimenMaterial describes the material obtained. If the LOINC test code also implicitly describes a material, this element may not conflict with the description. If desired, this component can provide a more detailed description of the material: LOINC codes only contain the materials at a main level.  This is in line with the agreements made in the IHE/Nictiz program e-Lab.  If the test is carried out on derived material (such as plasma), this element will still contain the material drawn (in this case, blood). In this case, the LOINC code will generally refer to plasma.	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.16	
<b>DCM::DefinitionCode</b>	SNOMED CT: 370133003 Specimen substance	
<b>DCM::ExampleValue</b>	Urine	
<b>DCM::ValueSet</b>	SpecimenMaterialCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.6
<b>Opties</b>		

<b>«data»</b>	Microorganism	
<b>Definitie</b>	In particular in microbiological determinations the subject of the test is an isolate of certain microorganism rather than a material. This concept provides the ability to capture information about this microorganism.	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.22	
<b>DCM::ValueSet</b>	MicroorganismCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.10
<b>Opties</b>		

<b>«data»</b>	CollectedVolume
<b>Definitie</b>	Total volume of the collected material. If it is necessary to determine the absolute amount of a particular substance in the collected material, the volume thereof must be given.
<b>Datatype</b>	PQ
<b>DCM::ConceptId</b>	NL-CM:13.1.23
<b>Opties</b>	

<b>«data»</b>	CollectionPeriod
<b>Definitie</b>	If the material has not been collected at a single point in time but over a certain period, this period can be captured in this concept. An example is 24 hour urine.
<b>Datatype</b>	TS
<b>DCM::ConceptId</b>	NL-CM:13.1.24
<b>Opties</b>	

<b>«data»</b>	CollectionDateTime
<b>Definitie</b>	Time at which the material was collected.
<b>Datatype</b>	TS
<b>DCM::ConceptId</b>	NL-CM:13.1.17
<b>DCM::DefinitionCode</b>	SNOMED CT: 399445004 specimen collection date
<b>DCM::ExampleValue</b>	10-07-2012 17:20
<b>Opties</b>	

<b>«data»</b>	ReceivedDateTime
<b>Definitie</b>	Date and time that the material is handed over at the laboratory or specimen collection center. This is the issue with material that is collected by the patient himself.
<b>Datatype</b>	TS
<b>DCM::ConceptId</b>	NL-CM:13.1.25
<b>Opties</b>	

<b>«data»</b>	CollectionMethod
<b>Definitie</b>	If relevant for the results, the method of obtaining the specimen can be entered as well.
<b>Datatype</b>	CD
<b>DCM::ConceptId</b>	NL-CM:13.1.18
<b>DCM::DefinitionCode</b>	SNOMED CT: 118171006 Specimen procedure
<b>DCM::ExampleValue</b>	Midstream
<b>DCM::ValueSet</b>	CollectionMethodCodelist
	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.2
<b>Opties</b>	

<b>«data»</b>	AnatomicalLocation
<b>Definitie</b>	Anatomic location where the material is collected, e.g. elbow
<b>Datatype</b>	CD

<b>DCM::ConceptId</b>	NL-CM:13.1.26	
<b>DCM::DefinitionCode</b>	SNOMED CT: 405814001 Procedure site - indirect	
<b>DCM::ValueSet</b>	SpecimenAnatomicalLocatio nCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.11
<b>Opties</b>		

«data»	Laterality	
<b>Definitie</b>	Laterality adds information about body side to the anatomic location, e.g. left	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.27	
<b>DCM::DefinitionCode</b>	SNOMED CT: 272741003 Laterality	
<b>DCM::ValueSet</b>	LateralityCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.12
<b>Opties</b>		

«data»	Morphology	
<b>Definitie</b>	Morphology describes morphological abnormalities of the anatomical location where the material is taken, for example wound, ulcer.	
<b>Datatype</b>	CD	
<b>DCM::ConceptId</b>	NL-CM:13.1.28	
<b>DCM::DefinitionCode</b>	SNOMED CT: 118168003 Specimen source morphology	
<b>DCM::ValueSet</b>	MorphologyCodelist	OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.13
<b>Opties</b>		

«data»	SpecimenSource	
<b>Definitie</b>	If the material is not collected directly from the patient but comes from a patient-related object, e.g. a cathetertip, this source of material can be recorded here.	
<b>Datatype</b>	ST	
<b>DCM::ConceptId</b>	NL-CM:13.1.29	
<b>DCM::DefinitionCode</b>	SNOMED CT: 127454002 Device specimen	
<b>Opties</b>		

«data»	Comment	
<b>Definitie</b>	Comments on administering the test, such as drawing material after a (glucose) stimulus or taking medicine.	
<b>Datatype</b>	ST	
<b>DCM::ConceptId</b>	NL-CM:13.1.19	
<b>DCM::DefinitionCode</b>	LOINC: 48767-8 Annotation comment	
<b>DCM::ExampleValue</b>	Na (glucose)stimulus	
<b>Opties</b>		

<b>«document»</b>	SpecimenMaterialCodelist
-------------------	--------------------------

<b>Definie</b>		
<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.6	
<b>Opties</b>		
<b>MonstermateriaalCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.6</b>
Codes	Coding Syst. Name	Coding System OID
SNOMED CT: < 105590001  substantie	SNOMED CT	2.16.840.1.113883.6.96

Specifieker waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

<b>«document»</b>	ContainerTypeCodelist	
<b>Definitie</b>		
<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.9	
<b>Opties</b>		
<b>ContainerTypeCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.9</b>
Codes	Coding Syst. Name	Coding System OID
SNOMED CT: < 434711009  Specimen container (physical object)	SNOMED CT	2.16.840.1.113883.6.96

Specifieker waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

<b>«document»</b>	MicroorganismCodelist	
<b>Definitie</b>		
<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.10	
<b>Opties</b>		
<b>MicroorganismeCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.10</b>
Codes	Coding Syst. Name	Coding System OID
SNOMED CT: ^2581000146104  simple referentieset voor micro-organismen (foundation metadata concept)	SNOMED CT	2.16.840.1.113883.6.96

<b>«document»</b>	CollectionMethodCodelist	
<b>Definitie</b>		
<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.2	
<b>Opties</b>		
<b>AfnameprocedureCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.2</b>
Codes	Coding Syst. Name	Coding System OID

SNOMED CT: <7138002   verrichting (verrichting) (alle concepten onder procedure)	SNOMED CT	2.16.840.1.113883.6.96
----------------------------------------------------------------------------------------	-----------	------------------------

Specifieker waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

«document»	SpecimenAnatomicalLocationCodelist	
<b>Definitie</b>		
<b>Datatype</b>		
<b>DCM::ValueSetBinding</b>	Extensible	
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.11	
<b>Opties</b>		
<b>MonsterAnatomischeLocatieCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.11</b>
Codes	Coding Syst. Name	Coding System OID
SNOMED CT: < 442083009  Anatomical or acquired body structure	SNOMED CT	2.16.840.1.113883.6.96

«document»	LateralityCodelist						
<b>Definitie</b>							
<b>Datatype</b>							
<b>DCM::ValueSetBinding</b>	Extensible						
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.12						
<b>Opties</b>							
<b>LateraliteitCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.12</b>					
Concept Name	Concept Code	CodeSys. Name	CodeSystem OID	Description			
Left	7771000	SNOMED CT	2.16.840.1.113883.6.96	Links			
Right	24028007	SNOMED CT	2.16.840.1.113883.6.96	Rechts			
Right and left	51440002	SNOMED CT	2.16.840.1.113883.6.96	Rechts en links			

«document»	MorphologyCodelist				
<b>Definitie</b>					
<b>Datatype</b>					
<b>DCM::ValueSetBinding</b>	Extensible				
<b>DCM::ValueSetId</b>	2.16.840.1.113883.2.4.3.11. 60.40.2.13.1.13				
<b>Opties</b>					
<b>MorfologieCodelijst</b>		<b>OID: 2.16.840.1.113883.2.4.3.11.60.40.2.13.1.13</b>			
Codes	Coding Syst. Name	Coding System OID			
SNOMED CT: < 49755003  Morphologically abnormal structure	SNOMED CT	2.16.840.1.113883.6.96			

Specifieker waardenlijsten zullen op termijn opgeleverd worden door het project "Lab-Terminologie"

## 1.8 Example Instances

LaboratoriumUitslag										
Resultaat Type	Resultaat Status	Monster		LaboratoriumTest						
		Monster materiaal	Afname DatumTijd	TestCode	Test DatumTijd	TestUitslag	Referentie Ondergrens	Referentie Bovengrens	Interpretatie Vlaggen	
Klinische chemie	Definitief	Bloed	12-06-2012 09:00	Natrium	12-06-2012 13:15	138 mmol/l	136 mmol/l	146 mmol/l		

LaboratoriumUitslag										
Resultaat Type	Resultaat Status	Monster		LaboratoriumTest						
		Monster materiaal	Afname DatumTijd	TestCode	Test DatumTijd	TestUitslag	Referentie Ondergrens	Referentie Bovengrens	Interpretatie Vlaggen	
Klinische chemie	Definitief	Bloed	23-05-2012 08:08	Chloride	23-05-2012 12:00	109 mmol/l	99 mmol/l	108 mmol/l	Boven referentie-waarde	

LaboratoriumUitslag										
Resultaat Type	Resultaat Status	Monster		LaboratoriumTest						
		Monster materiaal	Afname DatumTijd	TestCode	Test DatumTijd	TestUitslag	Referentie Ondergrens	Referentie Bovengrens	Interpretatie Vlaggen	
Virologie	Definitief	Bloed	16-01-2012 08:00	Hepatitis A IgM	16-01-2012 10:12	Negatief				

## 1.9 Instructions

## 1.10 Interpretation

## 1.11 Care Process

## 1.12 Example of the Instrument

## 1.13 Constraints

## 1.14 Issues

## 1.15 References

1. Nederlandse Vereniging voor Medische Microbiologie (2010) *ELab en Evt*. [Online] Beschikbaar op: [http://www.nvmm.nl/ict/vereniging/werkgroepen\\_commissies/elab-en-evt](http://www.nvmm.nl/ict/vereniging/werkgroepen_commissies/elab-en-evt) [Geraadpleegd: 23 juli 2014].

## 1.16 Functional Model

## 1.17 Traceability to other Standards

## 1.18 Disclaimer

The Health and Care Information Models (a.k.a Clinical Building Block) has been made in collaboration with several different parties in healthcare. These parties asked Nictiz to manage good maintenance and

development of the information models. Hereafter, these parties and Nictiz are referred to as the collaborating parties. The collaborating parties paid utmost attention to the reliability and topicality of the data in these Health and Care Information Models. Omissions and inaccuracies may however occur. The collaborating parties are not liable for any damages resulting from omissions or inaccuracies in the information provided, nor are they liable for damages resulting from problems caused by or inherent to distributing information on the internet, such as malfunctions, interruptions, errors or delays in information or services provided by the parties to you or by you to the parties via a website or via e-mail, or any other digital means. The collaborating parties will also not accept liability for any damages resulting from the use of data, advice or ideas provided by or on behalf of the parties by means of the Health and Care Information Models. The parties will not accept any liability for the content of information in this Health and Care Information Model to which or from which a hyperlink is referred. In the event of contradictions in mentioned Health and Care Information Model documents and files, the most recent and highest version of the listed order in the revisions will indicate the priority of the documents in question. If information included in the digital version of a Health and Care Information Model is also distributed in writing, the written version will be leading in case of textual differences. This will apply if both have the same version number and date. A definitive version has priority over a draft version. A revised version has priority over previous versions.

## **1.19 Terms of Use**

The user may use the Health and Care Information Models without limitations. The copyright provisions in the paragraph concerned apply to copying, distributing and passing on the Health and Care Information Models.

## **1.20 Copyrights**

A Health and Care Information Model qualifies as a work within the meaning of Section 10 of the Copyright Act (Auteurswet). Copyrights protect the Health and Care Information Models and these rights are owned by the cooperating parties.

The user may copy, distribute and pass on the information in this Health and Care Information Model under the conditions that apply for Creative Commons license Attribution-NonCommercial-ShareAlike 3.0 Netherlands (CC BY-NCSA-3.0).

The content is available under Creative Commons Attribution-NonCommercial-ShareAlike 3.0 (see also <http://creativecommons.org/licenses/by-nc-sa/3.0/nl/>)

This does not apply to information from third parties that sometimes is used and / or referred to in a Health and Care Information Model, for example to an international medical terminology system. Any (copyright) rights that protect this information are not owned by the cooperating parties but by those third parties.